


# ESTElift®

P D O T H R E A D



*Key to the Beauty Aesthetic*



*ESTElift has various of products with World's First PDO Threads.*

*ESTElift is absorbable Polydioxanone (PDO) Thread  
that has High Quality and Safety Proven Products.*

*ESTElift enables Stronger Skin Lifting & Tightening,  
Whitening & Wrinkle Removal with Immediate Effects & Rapid Recovery.*

**ESTElift®**  
P D O   T H R E A D





*ESTElift works on the system of creating collagen and fibrotic reactions within the deep tissue which ultimately pull tissues back.*

*ESTElift thread lifting offers easy and safe way for bringing the youth and beauty back.*



# MONO

ESTELift Mono Thread is the basic type of PDO absorbable thread use for skin firming and collagen generation. It is recommended for first time users of thread lifting treatment and for combination with other applied threads.




NEEDLE TYPE : SHARP

| NEEDLE(mm) |        | PDO(mm) |           |
|------------|--------|---------|-----------|
| GAUGE      | LENGTH | LENGTH  | USP       |
|            | 15mm   | 15mm    |           |
|            | 25mm   | 30mm    |           |
|            | 38mm   | 50mm    |           |
| 29G        | 50mm   | 70mm    | 6-0       |
| 27G        | 60mm   | 90mm    | 5-0       |
|            | 90mm   | 140mm   |           |
|            | 12mm   | 15mm    |           |
| 30G / 31G  | 25mm   | 30mm    | 6-0 / 7-0 |
|            | 38mm   | 50mm    |           |

# SCREW

ESTELift Screw Thread is screw-shaped mono type of thread rolled up on a needle. The screw shape of thread is perfectly retained even after the injection into the skin. It has great effect on creating volume effect for wrinkles and dimpled areas of the face.



NEEDLE TYPE : SHARP

| NEEDLE(mm) |        | PDO(mm) |     |
|------------|--------|---------|-----|
| GAUGE      | LENGTH | LENGTH  | USP |
|            | 38mm   | 50mm    |     |
| 26G / 27G  | 50mm   | 70mm    | 5-0 |
|            | 60mm   | 90mm    |     |
|            | 90mm   | 140mm   |     |
| 30G        | 38mm   | 50mm    |     |
| 29G        | 50mm   | 70mm    | 6-0 |
|            | 60mm   | 90mm    |     |



# TWIN

ESTElift Twin has two strands of suture are twisted along the length of the needle twisted. Connective tissue fibers are filled up at the spaces between the sutures of the twisted strands, so that they can be more firmly attached to it. Twin applies for deep wrinkles (Forehead, wrinkles).



NEEDLE TYPE : SHARP

| NEEDLE(mm)  |       |        | PDO(mm) |     |
|---|-------|--------|---------|-----|
| TYPE  | GAUGE | LENGTH | LENGTH  | USP |
|  | 25G   | 38mm   | 50mm    | 6-0 |
|   | 26G   | 50mm   | 70mm    | 7-0 |
|   | 27G   | 60mm   | 90mm    |     |
|   |       | 90mm   | 140mm   |     |

# FC(Fine Cut)NEEDLE



Estelift FC(Fine Cut)Needle breaks the adhesion between the dermal layer and the epidermal layer which separates the layers. Therefore, the suture thread is easily accessible to the SMAS layer and the tension of the thread can be stronger. Fine Cut Needle's tip of the needle is cut horizontally and the blade on the back of the needle is precisely processed.

NEEDLE TYPE : FC TYPE

| NEEDLE(mm)  |       |        |        |      | PDO(mm)   |           |        |        |     |
|---|-------|--------|--------|------|---|-----------|--------|--------|-----|
| TYPE  | GAUGE | LENGTH | LENGTH | USP  | TYPE  | GAUGE     | LENGTH | LENGTH | USP |
|  | 23G   | 38mm   | 50mm   | 2-0  |  | 26G / 27G | 38mm   | 50mm   | 5-0 |
|   |       | 50mm   | 70mm   |      |   |           | 50mm   | 70mm   |     |
|   |       | 70mm   | 100mm  |      |   |           | 70mm   | 100mm  |     |
|   |       | 90mm   | 140mm  |      |   |           | 90mm   | 140mm  |     |
|   |       | 100mm  | 150mm  |      |   |           | 90mm   | 140mm  |     |
|  | 25G   | 38mm   | 50mm   | 4-0  |  | 30G       | 38mm   | 50mm   | 6-0 |
|   |       | 50mm   | 70mm   | 50mm |   |           | 70mm   |        |     |
|   |       | 60mm   | 90mm   | 60mm |   |           | 90mm   |        |     |
|  | 27G   | 90mm   | 140mm  | 5-0  |  | 29G       | 60mm   | 90mm   | 6-0 |
|   |       | 25mm   | 30mm   |      |   |           | 60mm   | 90mm   |     |
|   |       | 38mm   | 50mm   |      |   |           | 60mm   | 90mm   |     |
|   |       | 50mm   | 70mm   |      |   |           | 60mm   | 90mm   |     |
|  | 29G   | 38mm   | 50mm   | 6-0  |  | 25G       | 38mm   | 50mm   | 6-0 |
|   |       | 50mm   | 70mm   |      |   |           | 50mm   | 70mm   |     |
|   |       | 70mm   | 100mm  |      |   |           | 70mm   | 100mm  |     |
|   |       | 60mm   | 90mm   |      |   |           | 70mm   | 100mm  |     |
|   |       | 60mm   | 90mm   |      |   | 27G       | 90mm   | 140mm  |     |





### FC Needle

## PALM TREE PDO

Estelift Palm tree PDO is connected to the main thread on the inside of the needle and to the auxiliary thread on the outside of needle. Palm tree PDO thread has one strand of thick suture that makes a large loop shape, thin strand of sutures are dispersed in the loop.



NEEDLE TYPE : FC TYPE

| NEEDLE(mm) |        | STEM PDO(mm) |     |          | BRANCH PDO(mm) |     |          |
|------------|--------|--------------|-----|----------|----------------|-----|----------|
| GAUGE      | LENGTH | LENGTH       | USP | Quantity | LENGTH         | USP | Quantity |
| 21G        | 38mm   | 50mm         | 4-0 | 1        | 50mm           | 5-0 | 14       |
|            |        | 175mm        |     |          |                |     |          |
|            | 60mm   | 80mm         |     |          |                |     |          |
|            |        | 225mm        |     |          |                |     |          |

## VINE SPRING

Estelift Vine spring thread is closely spaced dense coils that has the dense tunnel holding coil shaped sutures. It is upgraded screw-type PDO thread that gives natural volumizing effect for more dramatic result. Surgical site: Fallen cheeks and nasolabial fold.



NEEDLE TYPE : L TYPE

| NEEDLE(mm) |        | PDO(mm) |     |
|------------|--------|---------|-----|
| GAUGE      | LENGTH | LENGTH  | USP |
| 21G        | 60mm   | 340mm   | 4-0 |
| 23G        |        |         |     |
| 25G        |        | 300mm   |     |
| 27G        | 50mm   | 220mm   | 5-0 |
| 29G        |        | 6-0     |     |





## EYE CANNULA

Estelift Eye PDO thread has USP 6-0 which is mainly chosen in EYE / C, L for reducing the sensitive pain and swelling around eyes while treatment, It is inserted in the dermis layer under the skin around the eyes to improve elasticity and fine wrinkles of eye area.



NEEDLE TYPE : C / L TYPE

| SHAPE         | NEEDLE (mm) |        | PDO (mm) |     | USP |
|---------------|-------------|--------|----------|-----|-----|
|               | GAUGE       | LENGTH | LENGTH   |     |     |
| <i>C TYPE</i> | 30G         | 25mm   | 30mm     | 6-0 |     |
|               |             | 38mm   | 50mm     |     |     |
| <i>L TYPE</i> | 30G         | 25mm   | 30mm     | 6-0 |     |
|               |             | 38mm   | 50mm     |     |     |

## NOSE CANNULA

Estelift Nose thread is designed to optimize the nose tip and nose to correct. Cog thread is injected into dermal layer of the nasal septum. Suture thread is injected into the septum of the nose and is absorbed and hydrolyzed for long period of time. It enables maintaining perfect shape of the nose.



NEEDLE TYPE : N / L TYPE

| SHAPE         | NEEDLE (mm) |        | PDO (mm) |     |         |
|---------------|-------------|--------|----------|-----|---------|
|               | GAUGE       | LENGTH | LENGTH   | USP | PATTERN |
| <i>N TYPE</i> | 19G         | 38mm   | 30mm     | 2   |         |
|               |             |        | 40mm     |     |         |
|               |             |        | 50mm     |     |         |
|               |             |        | 60mm     |     |         |
| <i>L TYPE</i> | 19G         | 50mm   | 70mm     | 3-0 |         |
|               |             |        | 70mm     |     |         |
|               |             |        | 80mm     |     |         |
| <i>L TYPE</i> | 23G         | 38mm   | 50mm     | 3-0 |         |
|               |             | 50mm   | 70mm     |     |         |





## COG SHARP

Estelift COG Sharp is used for suture of tissue, ligation and fix by way of COG shape projecting. It can be injected without surgical incision. It has one way or two way (bi-directional) of COG projections. It usually used to surgical procedures correcting the jaw line.



NEEDLE TYPE : SHARP 

| NEEDLE(mm) |        | PDO(mm) |     |   |
|------------|--------|---------|-----|---|
| GAUGE      | LENGTH | LENGTH  | USP | PATTERN   |
| 23G        | 50mm   | 70mm    | 3-0 |  |
|            | 70mm   | 90mm    |     |   |
|            | 70mm   | 110mm   |     |   |
|            | 90mm   | 150mm   |     |  |
|            | 90mm   | 170mm   |     |   |
|            |        |         |     |   |

## COG SHARP 2D

Estelift COG Sharp is used for suture of tissue, ligation and fix by way of COG shape projecting. It can be injected without surgical incision. It has one way or two way (bi-directional) of COG projections. It usually used to surgical procedures correcting the jaw line.



NEEDLE TYPE : SHARP 

| NEEDLE(mm) |        | PDO(mm) |     |  |
|------------|--------|---------|-----|--|
| GAUGE      | LENGTH | LENGTH  | USP | PATTERN  |
| 23G        | 50mm   | 70mm    | 3-0 |  |
|            | 70mm   | 90mm    |     |  |
|            | 70mm   | 110mm   |     |  |
|            | 70mm   | 150mm   |     |  |
|            | 90mm   | 170mm   |     |  |
|            |        |         |     |  |





# COG 2D

Estelift 2D is composed with 19G 100mm needle (W Type needle, L Type needle, N Type needle) that can be injected without surgical incision. It is easier to use than normal needle, because it reduces pain and swelling during the procedure. It is usually used to correcting the jaw line.




NEEDLE TYPE : W / L / N TYPE

| SHAPE  | NEEDLE (mm) |        | PDO (mm)       |     |                      |
|--|-------------|--------|----------------|-----|----------------------|
|  | GAUGE       | LENGTH | LENGTH         | USP | PATTERN              |
| <br>W TYPE CANNULA        | 19G         | 100mm  | 160mm          | 0   |                      |
| <br>L TYPE CANNULA        | 19G         | 100mm  | 160mm          | 0   | ←←←←←→→→→→<br>2D (D) |
| <br>N TYPE (Blunt Needle) | 19G         | 100mm  | 150mm<br>200mm | 1-0 |                      |

# COG 360°/6D

Estelift 360° 6D is composed with 19G 100mm needle (W Type needle, L Type needle, N Type needle) that can be injected without surgical incision. It is easier to use than normal needle, because it reduces pain and swelling during the procedure. It is usually used to correcting the jaw line.



NEEDLE TYPE : W / L / N TYPE


| SHAPE  | NEEDLE (mm) |        | PDO (mm) |       |                         |
|--|-------------|--------|----------|-------|-------------------------|
|  | GAUGE       | LENGTH | LENGTH   | USP   | PATTERN                 |
| <br>W TYPE CANNULA        | 19G         | 18G    | 100mm    | 160mm | 2-0                     |
|  |             | 60mm   | 110mm    | 1-0   |                         |
|  |             | 100mm  | 150mm    | 1-0   |                         |
| <br>L TYPE CANNULA        | 19G         | 100mm  | 200mm    | 1-0   | ←←←←←→→→→→<br>360° (6D) |
| <br>N TYPE (Blunt Needle) | 19G         | 100mm  | 150mm    | 1-0   |                         |





## SCREW COG

Estelift Screw COG has two advantages that patients can experience strengths of COG and Screw. The suture is extended outward from the hub so that the suture can be pulled out during the procedure. It can prevent suture's rolling down and is used for all-round facial correction.




NEEDLE TYPE : SHARP

| NEEDLE(mm) |        | PDO(mm) |     |         |
|------------|--------|---------|-----|---------|
| GAUGE      | LENGTH | LENGTH  | USP | PATTERN |
| 26G        | 38mm   | 150mm   | 29G |         |
|            | 60mm   | 190mm   |     |         |
| 27G        | 38mm   | 70mm    | 5-0 |         |
|            | 38mm   | 80mm    |     |         |
|            | 60mm   | 110mm   |     |         |
|            | 90mm   | 140mm   |     |         |

## PRESSED COG

Estelift Pressed COG is new concept product, has the strongest lifting effect on the procedure. Pressed COG is shaped in compression type by shape processing (Squeezed suture). Tensile strength is more than three times compare to cutting cog. It prevents the breakage of the thread after the procedure by the facial muscle movement. It is mainly used for correcting the jaw line. It activates fibroblast and help the skin elasticity recovery and is slowly hydrolyzed during 6 to 8 months in the body.



REMARKS

Needle Option : N / L / W Type

| NEEDLE(mm)    |       |        | PDO(mm) |     |         |
|---------------|-------|--------|---------|-----|---------|
| SHAPE         | GAUGE | LENGTH | LENGTH  | USP | PATTERN |
|               | 19G   | 100mm  | 200mm   | 2   |         |
| <i>N TYPE</i> |       |        |         |     |         |
| <i>L TYPE</i> |       |        |         |     |         |
| <i>W TYPE</i> |       |        |         |     |         |



| TYPE                  | PACKING UNIT |        | STORAGE                     |
|-----------------------|--------------|--------|-----------------------------|
|                       | 1PACK        | 1 BOX  |                             |
| <b>Mono</b>           | 10PCS        | 100PCS | Room Temperature<br>(1~30°) |
| <b>Screw</b>          | 10PCS        | 100PCS |                             |
| <b>Twin</b>           | 10PCS        | 100PCS |                             |
| <b>Palm tree</b>      | 4PCS         | 40PCS  |                             |
| <b>Coil Spring</b>    | 2PCS         | 20PCS  |                             |
| <b>Mono Eye (C,L)</b> | 10PCS        | 100PCS |                             |
| <b>Nose COG</b>       | 4PCS         | 40PCS  |                             |
| <b>COG</b>            | 4PCS         | 40PCS  |                             |
| <b>Pressed COG</b>    | 4PCS         | 40PCS  |                             |
| <b>Screw COG</b>      | 4PCS         | 40PCS  |                             |

**ESTElift®**  
PDO THREAD



[www.ewellmedi.com](http://www.ewellmedi.com)

[www.estelift.com](http://www.estelift.com)

E-mail : [info@estelift.com](mailto:info@estelift.com)

Tel : +82 70 5089 3335