

# ESTElift®

P D O T H R E A D

*Key to the Beauty Aesthetic*





*ESTElift works on the system of creating collagen and fibrotic reactions within the deep tissue which ultimately pull tissues back.*

*ESTElift thread lifting offers easy and safe way for bringing the youth and beauty back.*



*ESTELIFT has various of products with World's First PDO Threads.*

*ESTELIFT is absorbable Polydioxanone (PDO) Thread that has High Quality and Safety Proven Products.*

*ESTELIFT enables Stronger Skin Lifting & Tightening, Whitening & Wrinkle Removal with Immediate Effects & Rapid Recovery.*



## MONO

ESTELIFT Mono Thread is the basic type of PDO absorbable thread use for skin firming and collagen generation. It is recommended for first time users of thread lifting treatment and for combination with other applied threads.



NEEDLE TYPE : SHARP

NEEDLE(mm)		PDO(mm)		
GAUGE	LENGTH	LENGTH	USP	
29G	15mm	15mm		6-0
	25mm	30mm		
	38mm	50mm		
	50mm	70mm		
	60mm	90mm	5-0	
	90mm	140mm		
30G / 31G	12mm	15mm		6-0 / 7-0
	25mm	30mm		
	38mm	50mm		



## SCREW

ESTElift Screw Thread is screw-shaped mono type of thread rolled on a needle. The screw shape of thread is perfectly retained even after the injection into the skin. It has great effect on creating volume effect for wrinkles and dimpled areas of the face.



NEEDLE TYPE : SHARP

GAUGE	NEEDLE(mm)		USP
	LENGTH	PDO(mm)	
26G / 27G	38mm	50mm	5-0
	50mm	70mm	
	60mm	90mm	
	90mm	140mm	
29G	38mm	50mm	6-0
	50mm	70mm	
	60mm	90mm	

## DOUBLE SCREW

ESTElift Double Screw thread is a screw-shaped multi-thread. The effect is twice stronger than injecting screw thread. The shape of the thread is perfectly retained even after the injection into the skin. It has an enhanced effect on creating volume effect for wrinkles and dimpled area of the face.



NEEDLE TYPE : SHARP

GAUGE	NEEDLE(mm)		USP
	LENGTH	PDO(mm)	
25 G	38mm	50mm	6-0
26 G	50mm	70mm	
	60mm	90mm	
27 G	90mm	140mm	7-0
	29G	38mm	
50mm		70mm	
60mm		90mm	





# TWIN

ESTELIFT Twin has two strands of suture are twisted along the length of the needle twisted. Connective tissue fibers are filled up at the spaces between the sutures of the twisted strands, so that they can be more firmly attached to it. Twin applies for deep wrinkles (Forehead, wrinkles).



NEEDLE TYPE : SHARP

NEEDLE(mm)			PDO(mm)	
TYPE	GAUGE	LENGTH	LENGTH	USP
	25G	38mm	50mm	6-0
	26G	29G	50mm	70mm
	27G	60mm	90mm	7-0
<b>TWIN</b>		90mm	140mm	

# FC(Fine Cut)NEEDLE



Estelift FC(Fine Cut)Needle breaks the adhesion between the dermal layer and the epidermal layer which separates the layers. Therefore, the suture thread is easily accessible to the SMAS layer and the tension of the thread can be stronger. Fine Cut Needle's tip of the needle is cut horizontally and the blade on the back of the needle is precisely processed.

NEEDLE TYPE : FC TYPE

NEEDLE(mm)					PDO(mm)				
TYPE	GAUGE	LENGTH	LENGTH	USP	TYPE	GAUGE	LENGTH	LENGTH	USP
	23G	38mm	50mm	2-0		26G / 27G	38mm	50mm	5-0
		50mm	70mm				50mm	70mm	
		70mm	100mm				70mm	100mm	
		90mm	140mm				90mm	140mm	
		100mm	150mm						
<b>FC MONO</b>	25G	38mm	50mm	4-0		25G / 26G / 27G	38mm	50mm	6-0
	26G	50mm	70mm	5-0			50mm	70mm	
	27G	60mm	90mm	5-0			60mm	90mm	
29G	30G	25mm	30mm	6-0		25G / 26G / 27G	38mm	50mm	6-0
		38mm	50mm				50mm	70mm	
		50mm	70mm				70mm	100mm	
		60mm	90mm				90mm	140mm	





### FC Needle

## PALM TREE PDO

Estelift Palm tree PDO is connected to the main thread on the inside of the needle and to the auxiliary thread on the outside of needle. Palm tree PDO thread has one strand of thick suture that makes a large loop shape, thin strand of sutures are dispersed in the loop.



NEEDLE TYPE : FC TYPE 

NEEDLE(mm)		STEM PDO(mm)			BRANCH PDO(mm)		
GAUGE	LENGTH	LENGTH	USP	Quantity	LENGTH	USP	Quantity
21G	38mm	50mm	4-0	1	50mm	5-0	14
		175mm					
	60mm	80mm					
		225mm					

## VINE SPRING

Estelift Vine spring thread is closely spaced dense coils that has the dense tunnel holding coil shaped sutures. It is upgraded screw-type PDO thread that gives natural volumizing effect for more dramatic result. Surgical site: Fallen cheeks and nasolabial fold.



NEEDLE TYPE : L TYPE 

NEEDLE(mm)		PDO(mm)	
GAUGE	LENGTH	LENGTH	USP
21G	60mm	340mm	4-0
23G		300mm	5-0
25G			
27G	50mm	220mm	6-0
29G			





## EYE CANNULA

Estelift Eye PDO thread has USP 6-0 which is mainly chosen in EYE / C, L for reducing the sensitive pain and swelling around eyes while treatment, It is inserted in the dermis layer under the skin around the eyes to improve elasticity and fine wrinkles of eye area.



NEEDLE TYPE : C / L TYPE

SHAPE	NEEDLE(mm)		PDO(mm)		USP
	GAUGE	LENGTH	LENGTH		
<i>C TYPE</i>	30G	25mm	30mm	6-0	
		38mm	50mm		
<i>L TYPE</i>	30G	25mm	30mm	6-0	
		38mm	50mm		

## NOSE CANNULA

Estelift Nose thread is designed to optimize the nose tip and nose to correct. Cog thread is injected into dermal layer of the nasal septum. Suture thread is injected into the septum of the nose and is absorbed and hydrolyzed for long period of time. It enables maintaining perfect shape of the nose.



NEEDLE TYPE : N / L TYPE

SHAPE	NEEDLE(mm)		PDO(mm)		
	GAUGE	LENGTH	LENGTH	USP	PATTERN
<i>N TYPE</i>	19G	38mm	30mm	2	
			40mm		
			50mm		
<i>L TYPE</i>	19G	50mm	60mm	3-0	
			70mm		
			80mm		
<i>L TYPE</i>	23G	38mm	50mm	3-0	
		50mm	70mm		





## COG SHARP

Estelift COG Sharp is used for suture of tissue, ligation and fix by way of COG shape projecting. It can be injected without surgical incision. It has one way or two way (bi-directional) of COG projections. It usually used to surgical procedures correcting the jaw line.



NEEDLE TYPE : SHARP 

NEEDLE(mm)		PDO(mm)		
GAUGE	LENGTH	LENGTH	USP	PATTERN
23G	50mm	70mm	3-0	
	70mm	90mm		
	70mm	110mm		
	90mm	150mm		
	90mm	170mm		

## COG SHARP 2D

Estelift COG Sharp is used for suture of tissue, ligation and fix by way of COG shape projecting. It can be injected without surgical incision. It has one way or two way (bi-directional) of COG projections. It usually used to surgical procedures correcting the jaw line.



NEEDLE TYPE : SHARP 

NEEDLE(mm)		PDO(mm)		
GAUGE	LENGTH	LENGTH	USP	PATTERN
23G	50mm	70mm	3-0	
	70mm	90mm		
	70mm	110mm		
	90mm	150mm		
	90mm	170mm		





# COG 2D

Estelift 2D is composed with 19G 100mm needle (W Type needle, L Type needle, N Type needle) that can be injected without surgical incision. It is easier to use than normal needle, because it reduces pain and swelling during the procedure. It is usually used to correcting the jaw line.



NEEDLE TYPE : W / L / N TYPE

SHAPE	NEEDLE(mm)		PDO(mm)		
	GAUGE	LENGTH	LENGTH	USP	PATTERN
 W TYPE CANNULA	19G	100mm	160mm	0	
 L TYPE CANNULA	19G	100mm	160mm	0	←←←←←→→→→→ 2D (D)
 N TYPE (Blunt Needle)	19G	100mm	150mm 200mm	1-0	

# COG 360°/6D

Estelift 360° 6D is composed with 19G 100mm needle (W Type needle, L Type needle, N Type needle) that can be injected without surgical incision. It is easier to use than normal needle, because it reduces pain and swelling during the procedure. It is usually used to correcting the jaw line.



NEEDLE TYPE : W / L / N TYPE


SHAPE	NEEDLE(mm)		PDO(mm)		
	GAUGE	LENGTH	LENGTH	USP	PATTERN
 W TYPE CANNULA	19G	100mm	160mm	2-0	
		60mm	110mm	1-0	
		100mm	150mm		
 L TYPE CANNULA	19G	100mm	200mm	1-0	←←←←←→→→→→ 360° (6D)
 N TYPE (Blunt Needle)	19G	100mm	150mm	1-0	





## SCREW COG

Estelift Screw COG has two advantages that patients can experience strengths of COG and Screw. The suture is extended outward from the hub so that the suture can be pulled out during the procedure. It can prevent suture's rolling down and is used for all-round facial correction.



NEEDLE TYPE : SHARP

NEEDLE(mm)		PDO(mm)		
GAUGE	LENGTH	LENGTH	USP	PATTERN
26G	38mm	150mm	29G	
	60mm	190mm		
27G	38mm	70mm	5-0	
	38mm	80mm		
	60mm	110mm		
	90mm	140mm		

## MOLDING COG

Estelift Pressed COG is new concept product, has the strongest lifting effect on the procedure. Pressed COG is shaped in compression type by shape processing (Squeezed suture). Tensile strength is more than three times compare to cutting cog. It prevents the breakage of the thread after the procedure by the facial muscle movement.


It is mainly used for correcting the jaw line. It activates fibroblast and help the skin elasticity recovery and is slowly hydrolyzed during 6 to 8 months in the body.

Needle Option:




NEEDLE(mm)		PDO(mm)	
GAUGE	LENGTH	LENGTH	USP
18G	60mm ~ 100mm	150mm ~ 200mm	2
19G			







*PCJL / PCJW*




*PCKL / PCKW*




*PCHL / PCHW*




*PCIL / PCIW*




*PCML / PCMW*



*PCNL / PCNW*



*PCLL / PCLW*



*PCFL / PCFW*



TYPE	PACKING UNIT		STORAGE
	1 PACK	1 BOX	
<b>Mono</b>	10PCS	100PCS	Room Temperature (1~30°)
<b>Screw / Double Screw</b>	10PCS	100PCS	
<b>Twin</b>	10PCS	100PCS	
<b>Palm tree</b>	4PCS	40PCS	
<b>Vine Spring</b>	2PCS	20PCS	
<b>Mono Eye (C,L)</b>	10PCS	100PCS	
<b>Nose COG</b>	4PCS	40PCS	
<b>COG</b>	4PCS	40PCS	
<b>Molding COG</b>	4PCS	40PCS	
<b>Screw COG</b>	4PCS	40PCS	

**ESTElift®**  
P O O T H R E A D



[www.estelift.com](http://www.estelift.com)  
E-mail : [info@estelift.com](mailto:info@estelift.com)  
Tel : +82 70 5089 3335